

Printer	Language	InBody270 - Firmware needed	InBody370s - Firmware needed	InBody570 - Firmware needed	InBody770 - Firmware needed	InBody970 - Firmware needed	InBodyS10 DNO1 - Firmware need	InBodyS10 - Firmware needed	InBody8WA2.0 - Firmware needed
Samsung M2026W	SPL2009-600	O	O	O	O	O	210	O	O
Brother HL1114	Brother GDI	422 (Kernel 1006)	422 (Kernel 1006)	422 (Kernel 1006)	422 (Kernel 1006)	426 (Kernel K003)	X	X - Requested	426 (Kernel K003)
Brother HL1210W*	Brother GDI	422 (Kernel 1006)	422 (Kernel 1006)	422 (Kernel 1006)	422 (Kernel 1006)	422 (Kernel K003)	X	X - Requested	422 (Kernel K003)
Brother HL-2270DW	PCL	O	O	O	O	O	Unsure when added	O	O
Brother HL-5100DN	PCL	O	O	O	O	O	Unsure when added	O	O
Brother HL-L2370D	PCL	O	O	O	O	O	Unsure when added	O	O
Brother HL-L3210CW	PCL	O	O	O	O	O	Unsure when added	O	O
Brother HL-L3230CDW	PCL	O	O	O	O	O	Unsure when added	O	O
HP - Envy5000	PCL	O	O	O	O	O	Unsure when added	O	O
HP ENVY 6075	PCL3GUI	382	381	381	381	380	X	390	381
HP ENVY Pro 6475	PCL3GUI	382	381	381	381	380	X	390	381
HP Laser 107A*	SPL2011	O	O	O	O	O	201	O	O
HP LaserJet M110w*	PCLms	400	400	400	400	405	X	390	405
HP LaserJet M110WE*	PCLms	400	400	400	400	405	X	390	405
HP LaserJet Pro M102w	PCLms	400	400	400	400	405	X	390	405
HP LaserJet Pro M15w	PCLms	400	400	400	400	405	X	390	405
HP LaserJet Pro M203dw	PCLms	400	400	400	400	405	X	390	405
HP LaserJet Pro M209dw	PCLms	400	400	400	400	405	X	390	405
HP LaserJet Pro M29w	PCLms	400	400	400	400	405	X	390	405
HP Neverstop Laser 1001nw	PCLms	400	400	400	400	405	X	390	405
HP OfficeJet 6950	PCL	O	O	O	O	O	Unsure when added	O	O
HP OfficeJet Pro 6230*	PCL3GUI	382	381	381	381	380	X	390	381
Pantum P3300 Series	PCL	O	O	O	O	O	Unsure when added	O	O
Samsung M2020*	SPL2009-600	O	O	O	O	O	210	O	O
Samsung ML-2165W	SPL2009-600	O	O	O	O	O	210	O	O
Samsung ML-2955ND	SPL2009-600	O	O	O	O	O	210	O	O

O - Available any firmware
X - Unavailable
* tested by InBody Europe