

Date: July 31st, 2023

Letter of Notification

To whom it may concern,

InBody Co., Ltd. hereby certifies that:

For **970DM-0430 Digital Firmware / 970BWA_MS-A022**,

Does not require any Loadcell or Impedance Calibration after firmware upgrade.
(*above firmware version **DNO3-970DM-0383**)

Note: Due to the potential for defects to occur during firmware upgrade, there is a possibility of requiring replaceable Parts or Calibration.

Confirmed by
Sam Yoo
CS Team Leader
InBody Co., Ltd.