

Document Number : IB-24-17

Released Date : 2024.07.15

Recipient : Domestic sales, Overseas sales

Reference : Quality Part, Production Part


Title : Removal of the Tapes for LCD FFC Cable and LNO in InBody270

1. We are much appreciated for the hard work from branches and subsidiaries.

2. Change Classification (Reason)

The Insulation Tape in the back of the LNO Board and the EMI tape of the LCD FFC Cable have been removed due to the improvement of LNO6 Board EMC.

3. Changed Details

Before	After
	
	

4. Date of Implementation

Model	S/N	Manufacture Date
InBody270	F924003417	2024-06-10

Headquarter : InBody Bldg., 625, Eonju-ro, Gangnam-gu, Seoul 06106 Korea

Manufacturer : 15, Heugam-gil, Ipjang-myeon Seobuk-gu, Cheonan-si Chungcheongnam-do 31025 KOREA

Tel : 010-9909-2372

Fax : (041)581-3103

Department in charge : Quality Management Team of Quality Part

Person in charge : TaeGu Kang E-mail : kp016250@inbody.com

IB-105-035

InBody Co., Ltd

Rev.0