

Document Number : IB-24-19

Released Date : 2024.07.22

Recipient : Domestic sales, Overseas sales

Reference : Quality Part, Production Part

Title : Change in InBody970 Digital Program (0460ver)

1. We are much appreciated for the hard work from branches and subsidiaries.

2. Change Classification (Reason)

- Modifying UI and Function

3. Changed Details

- 1) Application of multilingual text for the face maker on the bottom of the main screen
- 2) Addition of outputs, 'Adjusted SMI' and 'Adjusted FFM' (A4, Thermal)
- 3) Terminology change in Taiwanese (Traditional Chinese)

4. Date of Implementation

Model	Program	S/N	Manufacture Date
InBody970	Digital Program(0460Ver)	I424000328	2024/7/1