

Document Number : IB-24-26

Released Date : 2024.09.12

Recipient : Domestic sales, Overseas sales

Reference : Quality Part, Production Part

Title : Change in InBodyJ30/J50 Motor Program (114ver)

1. We are much appreciated for the hard work from branches and subsidiaries.

2. Change Classification (Reason)

- Quality Improvement

3. Changed Details

- The motor program has been changed to address the issue that the previous measurement value was displayed when height measurement was completed using the pressure mode in the J30M-113ver. It has been corrected to ensure normal measurement even when using the pressure mode.

4. Date of Implementation

| Model     | Program     | S/N        | Application Date |
|-----------|-------------|------------|------------------|
| InBodyJ30 | J30M-114ver | C524000110 | 2024-08-27       |
| InBodyJ50 | J30M-114ver | C924000071 | 2024-08-27       |